

# NEMS-FTIR-TD WITH EMILIE™

## Nanoelectromechanical Fourier Transform Infrared Spectroscopy with in-situ Thermal Desorption

October 17, 2024

### 1 INTRODUCTION

Nanoelectromechanical Fourier Transform Infrared Spectroscopy (NEMS-FTIR) with EMILIE™ is a unique and powerful analytical tool which allows for in-situ thermal desorption (TD) analysis [1]. Unlike standard thermogravimetric method, the analyte is placed directly on the EMILIE™ sampling and sensing chip and monitored in-situ as the sensor is heated, allowing for thermal desorption of semi-volatile compounds and separation of mixtures of analytes [1].

### 2 METHODOLOGY

Theobromine (Sigma-Aldrich C0750-5G) was sampled from a liquid solution (10 µg/mL in distilled deionized water) using standard aerosol techniques as described in the EMILIE™ sampling guide. The aerosol was generated using a pneumatic jet nebulizer (Meinhard, CytoNeb) with a self-aspiration capillary (20 µL min<sup>-1</sup>) and pressurized dry air at 3 bar. A spray chamber (Meinhard, CytoSpray) was used for removal of larger droplets. The analyte flow was directed through a diffusion dryer (TOPAS DDU 570/L) prior to sampling by impaction on the EMILIE™ nanomechanical sampling and sensing chip using the EMILIE™ sampling accessory.

Measurements with EMILIE™ are performed at a pressure of  $p < 1 \times 10^{-3}$  mbar. The EMILIE™ nanomechanical sampling and sensing chip was subjected to twenty successive ten-minute

heating cycles at 80°C. After 10 minutes heating, a spectrum was recorded at 10°C and a new heating cycle was started. The procedure was reported until the theobromine sample was completely desorbed from the surface of the chip.

### 3 RESULTS




Figure 1: NEMS-FTIR spectrum of 43ng theobromine recorded with EMILIE™ .

Due to theobromine's semivolatility in vacuum, it undergoes sublimation and desorbs from the EMILIE™ nanomechanical sampling and sensing chip during the heating cycle. EMILIE™ allows for the simultaneous monitoring of the process by a combination of photothermal IR spectroscopy (Figure 1), and thermogravimetric analysis (Figure 2).

